








AngioGATE[®]
Introducers


Kimal AngioGATE[®] arterial sheath introducers offer superb performance and quality. The range of configurations offer an effective and reliable solution for a variety of percutaneous vascular procedures.


PRODUCT	DESCRIPTION	CODE	BOX QUANTITY																																																																														
	<p>HV Introducers (No guidewire)</p> <table border="1"> <thead> <tr> <th>Fr</th> <th>Length</th> <th>Wire diameter</th> <th>Wire core/Tip</th> <th>Wire length</th> <th>Needle</th> </tr> </thead> <tbody> <tr><td>4</td><td>12cm</td><td>Tipped to .038"</td><td>N/A</td><td>N/A</td><td>N/A</td></tr> <tr><td>5</td><td>12cm</td><td>Tipped to .038"</td><td>N/A</td><td>N/A</td><td>N/A</td></tr> <tr><td>6</td><td>12cm</td><td>Tipped to .038"</td><td>N/A</td><td>N/A</td><td>N/A</td></tr> <tr><td>7</td><td>12cm</td><td>Tipped to .038"</td><td>N/A</td><td>N/A</td><td>N/A</td></tr> <tr><td>8</td><td>12cm</td><td>Tipped to .038"</td><td>N/A</td><td>N/A</td><td>N/A</td></tr> <tr><td>9</td><td>12cm</td><td>Tipped to .038"</td><td>N/A</td><td>N/A</td><td>N/A</td></tr> <tr><td>5</td><td>12cm</td><td>Tipped to .038"</td><td>N/A</td><td>N/A</td><td>N/A</td></tr> <tr><td>6</td><td>12cm</td><td>Tipped to .038"</td><td>N/A</td><td>N/A</td><td>N/A</td></tr> <tr><td>7</td><td>12cm</td><td>Tipped to .038"</td><td>N/A</td><td>N/A</td><td>N/A</td></tr> </tbody> </table>	Fr	Length	Wire diameter	Wire core/Tip	Wire length	Needle	4	12cm	Tipped to .038"	N/A	N/A	N/A	5	12cm	Tipped to .038"	N/A	N/A	N/A	6	12cm	Tipped to .038"	N/A	N/A	N/A	7	12cm	Tipped to .038"	N/A	N/A	N/A	8	12cm	Tipped to .038"	N/A	N/A	N/A	9	12cm	Tipped to .038"	N/A	N/A	N/A	5	12cm	Tipped to .038"	N/A	N/A	N/A	6	12cm	Tipped to .038"	N/A	N/A	N/A	7	12cm	Tipped to .038"	N/A	N/A	N/A																				
Fr	Length	Wire diameter	Wire core/Tip	Wire length	Needle																																																																												
4	12cm	Tipped to .038"	N/A	N/A	N/A																																																																												
5	12cm	Tipped to .038"	N/A	N/A	N/A																																																																												
6	12cm	Tipped to .038"	N/A	N/A	N/A																																																																												
7	12cm	Tipped to .038"	N/A	N/A	N/A																																																																												
8	12cm	Tipped to .038"	N/A	N/A	N/A																																																																												
9	12cm	Tipped to .038"	N/A	N/A	N/A																																																																												
5	12cm	Tipped to .038"	N/A	N/A	N/A																																																																												
6	12cm	Tipped to .038"	N/A	N/A	N/A																																																																												
7	12cm	Tipped to .038"	N/A	N/A	N/A																																																																												
	<p>HV Introducers</p> <table border="1"> <thead> <tr> <th>Fr</th> <th>Length</th> <th>Wire diameter</th> <th>Wire core/Tip</th> <th>Wire length</th> <th>Needle</th> </tr> </thead> <tbody> <tr><td>4</td><td>12cm</td><td>0.035"</td><td>Dbl flex/3mm "J"</td><td>50cm</td><td>N/A</td></tr> <tr><td>5</td><td>12cm</td><td>0.038"</td><td>Dbl flex/3mm "J"</td><td>50cm</td><td>N/A</td></tr> <tr><td>6</td><td>12cm</td><td>0.038"</td><td>Dbl flex/3mm "J"</td><td>50cm</td><td>N/A</td></tr> <tr><td>7</td><td>12cm</td><td>0.038"</td><td>Dbl flex/3mm "J"</td><td>50cm</td><td>N/A</td></tr> <tr><td>8</td><td>12cm</td><td>0.038"</td><td>Dbl flex/3mm "J"</td><td>50cm</td><td>N/A</td></tr> <tr><td>9</td><td>12cm</td><td>0.038"</td><td>Dbl flex/3mm "J"</td><td>50cm</td><td>N/A</td></tr> <tr><td>4</td><td>24cm</td><td>0.035"</td><td>Dbl flex/3mm "J"</td><td>80cm</td><td>N/A</td></tr> <tr><td>5</td><td>24cm</td><td>0.038"</td><td>Dbl flex/3mm "J"</td><td>80cm</td><td>N/A</td></tr> <tr><td>6</td><td>24cm</td><td>0.038"</td><td>Dbl flex/3mm "J"</td><td>80cm</td><td>N/A</td></tr> <tr><td>7</td><td>24cm</td><td>0.038"</td><td>Dbl flex/3mm "J"</td><td>80cm</td><td>N/A</td></tr> <tr><td>8</td><td>24cm</td><td>0.038"</td><td>Dbl flex/3mm "J"</td><td>80cm</td><td>N/A</td></tr> <tr><td>8</td><td>24cm</td><td>0.038"</td><td>Dbl flex/3mm "J"</td><td>80cm</td><td>N/A</td></tr> </tbody> </table>	Fr	Length	Wire diameter	Wire core/Tip	Wire length	Needle	4	12cm	0.035"	Dbl flex/3mm "J"	50cm	N/A	5	12cm	0.038"	Dbl flex/3mm "J"	50cm	N/A	6	12cm	0.038"	Dbl flex/3mm "J"	50cm	N/A	7	12cm	0.038"	Dbl flex/3mm "J"	50cm	N/A	8	12cm	0.038"	Dbl flex/3mm "J"	50cm	N/A	9	12cm	0.038"	Dbl flex/3mm "J"	50cm	N/A	4	24cm	0.035"	Dbl flex/3mm "J"	80cm	N/A	5	24cm	0.038"	Dbl flex/3mm "J"	80cm	N/A	6	24cm	0.038"	Dbl flex/3mm "J"	80cm	N/A	7	24cm	0.038"	Dbl flex/3mm "J"	80cm	N/A	8	24cm	0.038"	Dbl flex/3mm "J"	80cm	N/A	8	24cm	0.038"	Dbl flex/3mm "J"	80cm	N/A		
Fr	Length	Wire diameter	Wire core/Tip	Wire length	Needle																																																																												
4	12cm	0.035"	Dbl flex/3mm "J"	50cm	N/A																																																																												
5	12cm	0.038"	Dbl flex/3mm "J"	50cm	N/A																																																																												
6	12cm	0.038"	Dbl flex/3mm "J"	50cm	N/A																																																																												
7	12cm	0.038"	Dbl flex/3mm "J"	50cm	N/A																																																																												
8	12cm	0.038"	Dbl flex/3mm "J"	50cm	N/A																																																																												
9	12cm	0.038"	Dbl flex/3mm "J"	50cm	N/A																																																																												
4	24cm	0.035"	Dbl flex/3mm "J"	80cm	N/A																																																																												
5	24cm	0.038"	Dbl flex/3mm "J"	80cm	N/A																																																																												
6	24cm	0.038"	Dbl flex/3mm "J"	80cm	N/A																																																																												
7	24cm	0.038"	Dbl flex/3mm "J"	80cm	N/A																																																																												
8	24cm	0.038"	Dbl flex/3mm "J"	80cm	N/A																																																																												
8	24cm	0.038"	Dbl flex/3mm "J"	80cm	N/A																																																																												

PRODUCT	DESCRIPTION	CODE	BOX QUANTITY
	HV Introducers RO Band		
	Fr Length		
	5 12cm RO band introducer	K63/0512B	10
	6 12cm RO band introducer	K63/0612B	10
	7 12cm RO band introducer	K63/0712B	10
	5 5.5cm RO band introducer kit	K63/055MAB	10
	6 5.5cm RO band introducer kit	K63/065MAB	10
	Micro-access introducers		
Fr Length Wire diameter Wire core/Tip Wire length Needle			
4 5.5cm 0.018" Steel/Steel 40cm 21G 4cm non echo	K63/045MAP	10	
5 5.5cm 0.018" Steel/Steel 40cm 21G 4cm non echo	K63/055MAP	10	
6 5.5cm 0.018" Steel/Steel 40cm 21G 4cm non echo	K63/065MAP	10	
5 5.5cm 0.018" Steel/Steel 40cm 21G 7cm non echo	K63/055MAB	10	
6 5.5cm 0.018" Steel/Steel 40cm 21G 7cm non echo	K63/065MAB	10	

PRODUCT	DESCRIPTION	CODE	BOX QUANTITY			
	MicroIntroducers					
	<table border="1"> <thead> <tr> <th data-bbox="947 419 1211 443">Needle</th> <th data-bbox="1211 419 1581 443">Guidewire Core/Tip</th> <th data-bbox="1581 419 1749 443">Sheath/Dilator</th> </tr> </thead> </table>	Needle	Guidewire Core/Tip	Sheath/Dilator		
	Needle	Guidewire Core/Tip	Sheath/Dilator			
	21g x 4cm Echogenic	40cm x .018" Nitinol/Platinum	4fr x 10cm	K63/41040NP4E	10	
	21g x 7cm Echogenic	40cm x .018" Nitinol/Platinum	4fr x 10cm	K63/41040NP7E	10	
	21g x 7cm Echogenic	40cm x .018" Nitinol/Platinum	5fr x 10cm	K63/51040NP7E	10	
	21g x 7cm Echogenic	40cm x .018" Nitinol/Platinum with stiffened dilator	4fr x 10cm	K63/41040NP7E/SD	10	
	21g x 7cm Echogenic	40cm x .018" Nitinol/Platinum with stiffened dilator	5fr x 10cm	K63/51040NP7E/SD	10	
	21g x 7cm Echogenic	40cm x .018" Nitinol/Steel	4fr x 10cm	K63/41040NS7	10	
	21g x 7cm Echogenic	40cm x .018" Nitinol/Steel	5fr x 10cm	K63/51040NS7E	10	
	21g x 7cm Echogenic	60cm x .018" Nitinol/Steel	4fr x 10cm	K63/41060NS7E	10	
	21g x 7cm Echogenic	60cm x .018" Nitinol/Steel	5fr x 10cm	K63/51060NS7E	10	
	21g x 7cm Smooth	40cm x .018" Steel/Steel	4fr x 10cm	K63/41040SS7S	10	
	21g x 7cm Smooth	40cm x .018" Steel/Steel	5fr x 10cm	K63/51040SS7S	10	
	21g x 7cm Smooth	40cm x .018" Steel/Steel	4fr x 10cm	K63/41040SS7E	10	
	21g x 7cm Echogenic	40cm x .018" Steel/Steel	5fr x 10cm	K63/51040SS7E	10	
	21g x 7cm Smooth	60cm x .018" Steel/Steel	4fr x 10cm	K63/41060SS7E	10	
	21g x 7cm Smooth	60cm x .018" Steel/Steel	5fr x 10cm	K63/51060SS7E	10	
	21g x 7cm Smooth	60cm x .018" Nitinol/Platinum	4fr x 10cm	K63/41060NP7E	10	
	21g x 7cm Smooth	60cm x .018" Nitinol/Platinum	5fr x 10cm	K63/51060NP7E	10	

PRODUCT	DESCRIPTION	CODE	BOX QUANTITY																																																																														
	Transradial Introducers																																																																																
	<table border="1"> <thead> <tr> <th>Fr</th> <th>Length</th> <th>Wire diameter</th> <th>Wire core/Tip</th> <th>Wire length</th> <th>Needle</th> </tr> </thead> <tbody> <tr><td>5</td><td>12cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>50cm</td><td>21G 4cm non echo</td></tr> <tr><td>6</td><td>12cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>50cm</td><td>21G 4cm non echo</td></tr> <tr><td>7</td><td>12cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>50cm</td><td>21G 4cm non echo</td></tr> <tr><td>4</td><td>12cm</td><td>0.018"</td><td>Steel/Steel</td><td>50cm</td><td>21G 4cm non echo</td></tr> <tr><td>5</td><td>12cm</td><td>0.018"</td><td>Steel/Steel</td><td>50cm</td><td>21G 4cm non echo</td></tr> <tr><td>6</td><td>12cm</td><td>0.018"</td><td>Steel/Steel</td><td>50cm</td><td>21G 4cm non echo</td></tr> <tr><td>5</td><td>19cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>80cm</td><td>21G 4cm non echo</td></tr> <tr><td>6</td><td>19cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>80cm</td><td>21G 4cm non echo</td></tr> <tr><td>7</td><td>19cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>80cm</td><td>21G 4cm non echo</td></tr> <tr><td>5</td><td>24cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>80cm</td><td>21G 4cm non echo</td></tr> <tr><td>6</td><td>24cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>80cm</td><td>21G 4cm non echo</td></tr> <tr><td>7</td><td>24cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>80cm</td><td>21G 4cm non echo</td></tr> </tbody> </table>	Fr	Length	Wire diameter	Wire core/Tip	Wire length	Needle	5	12cm	0.018"	Nitinol/Platinum	50cm	21G 4cm non echo	6	12cm	0.018"	Nitinol/Platinum	50cm	21G 4cm non echo	7	12cm	0.018"	Nitinol/Platinum	50cm	21G 4cm non echo	4	12cm	0.018"	Steel/Steel	50cm	21G 4cm non echo	5	12cm	0.018"	Steel/Steel	50cm	21G 4cm non echo	6	12cm	0.018"	Steel/Steel	50cm	21G 4cm non echo	5	19cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo	6	19cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo	7	19cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo	5	24cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo	6	24cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo	7	24cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo		
	Fr	Length	Wire diameter	Wire core/Tip	Wire length	Needle																																																																											
	5	12cm	0.018"	Nitinol/Platinum	50cm	21G 4cm non echo																																																																											
	6	12cm	0.018"	Nitinol/Platinum	50cm	21G 4cm non echo																																																																											
	7	12cm	0.018"	Nitinol/Platinum	50cm	21G 4cm non echo																																																																											
	4	12cm	0.018"	Steel/Steel	50cm	21G 4cm non echo																																																																											
	5	12cm	0.018"	Steel/Steel	50cm	21G 4cm non echo																																																																											
	6	12cm	0.018"	Steel/Steel	50cm	21G 4cm non echo																																																																											
	5	19cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo																																																																											
	6	19cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo																																																																											
	7	19cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo																																																																											
	5	24cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo																																																																											
	6	24cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo																																																																											
	7	24cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo																																																																											
	Transradial Introducers with Hydrophilic Coating																																																																																
	<table border="1"> <thead> <tr> <th>Fr</th> <th>Length</th> <th>Wire diameter</th> <th>Wire core/Tip</th> <th>Wire length</th> <th>Needle</th> </tr> </thead> <tbody> <tr><td>5</td><td>12cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>50cm</td><td>21G 4cm non echo</td></tr> <tr><td>6</td><td>12cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>50cm</td><td>21G 4cm non echo</td></tr> <tr><td>7</td><td>12cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>50cm</td><td>21G 4cm non echo</td></tr> <tr><td>5</td><td>12cm</td><td>0.018"</td><td>Steel/Steel</td><td>50cm</td><td>21G 4cm non echo</td></tr> <tr><td>6</td><td>12cm</td><td>0.018"</td><td>Steel/Steel</td><td>50cm</td><td>21G 4cm non echo</td></tr> <tr><td>5</td><td>19cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>50cm</td><td>21G 4cm non echo</td></tr> <tr><td>6</td><td>19cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>50cm</td><td>21G 4cm non echo</td></tr> <tr><td>5</td><td>24cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>80cm</td><td>21G 4cm non echo</td></tr> <tr><td>6</td><td>24cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>80cm</td><td>21G 4cm non echo</td></tr> <tr><td>7</td><td>24cm</td><td>0.018"</td><td>Nitinol/Platinum</td><td>80cm</td><td>21G 4cm non echo</td></tr> </tbody> </table>	Fr	Length	Wire diameter	Wire core/Tip	Wire length	Needle	5	12cm	0.018"	Nitinol/Platinum	50cm	21G 4cm non echo	6	12cm	0.018"	Nitinol/Platinum	50cm	21G 4cm non echo	7	12cm	0.018"	Nitinol/Platinum	50cm	21G 4cm non echo	5	12cm	0.018"	Steel/Steel	50cm	21G 4cm non echo	6	12cm	0.018"	Steel/Steel	50cm	21G 4cm non echo	5	19cm	0.018"	Nitinol/Platinum	50cm	21G 4cm non echo	6	19cm	0.018"	Nitinol/Platinum	50cm	21G 4cm non echo	5	24cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo	6	24cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo	7	24cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo														
	Fr	Length	Wire diameter	Wire core/Tip	Wire length	Needle																																																																											
	5	12cm	0.018"	Nitinol/Platinum	50cm	21G 4cm non echo																																																																											
	6	12cm	0.018"	Nitinol/Platinum	50cm	21G 4cm non echo																																																																											
	7	12cm	0.018"	Nitinol/Platinum	50cm	21G 4cm non echo																																																																											
	5	12cm	0.018"	Steel/Steel	50cm	21G 4cm non echo																																																																											
	6	12cm	0.018"	Steel/Steel	50cm	21G 4cm non echo																																																																											
	5	19cm	0.018"	Nitinol/Platinum	50cm	21G 4cm non echo																																																																											
	6	19cm	0.018"	Nitinol/Platinum	50cm	21G 4cm non echo																																																																											
	5	24cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo																																																																											
6	24cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo																																																																												
7	24cm	0.018"	Nitinol/Platinum	80cm	21G 4cm non echo																																																																												

PRODUCT	DESCRIPTION	CODE	BOX QUANTITY
	Peelaway Introducers		
	Sheath Peelaway 5fr x20cm System Kit	K57/0520K	5
	Sheath Peelaway 6fr x20cm System Kit	K57/0620K	5
	Sheath Peelaway 7fr x 20cm	K57/0720	5
	Sheath Peelaway 7fr x 20cm System Kit	K57/0720K	5
	Sheath Peelaway 8fr x 20cm	K57/0820	5
	Sheath Peelaway 8fr x 20cm System Kit	K57/0820K	5
	Sheath Peelaway 9fr x 20cm	K57/0920	5
	Sheath Peelaway 9fr x 20cm System Kit	K57/0920K	5
	Sheath Peelaway 10 fr x 20cm System Kit	K57/1020K	5
	Sheath Peelaway 10.5 fr x 20cm System Kit	K57/10520K	5
	Sheath, Peel Away 14 fr, with Needle	K57/11/SE	5
	Sheath Peelaway 11fr x 20cm System Kit	K57/1120K	5
	Sheath Peelaway 12fr x 20cm System Kit	K57/1220K	5
Sheath Peelaway 12.5fr x 20cm System Kit	K57/12520K	5	
Sheath Peelaway 13fr x 20cm System Kit	K57/1320K	5	

PRODUCT	DESCRIPTION	CODE	BOX QUANTITY																																																
	Vessel Dilator																																																		
	<table border="1"> <thead> <tr> <th>Fr</th> <th>Length</th> <th>Wire Diameter</th> </tr> </thead> <tbody> <tr> <td>4</td> <td>20cm</td> <td>025"</td> </tr> <tr> <td>4</td> <td>20cm</td> <td>035"</td> </tr> <tr> <td>5</td> <td>20cm</td> <td>035</td> </tr> <tr> <td>6</td> <td>20cm</td> <td>035"</td> </tr> <tr> <td>6</td> <td>20cm</td> <td>038"</td> </tr> <tr> <td>7</td> <td>20cm</td> <td>038"</td> </tr> <tr> <td>8</td> <td>20cm</td> <td>038"</td> </tr> <tr> <td>9</td> <td>20cm</td> <td>038"</td> </tr> <tr> <td>10</td> <td>20cm</td> <td>038"</td> </tr> <tr> <td>11</td> <td>20cm</td> <td>038"</td> </tr> <tr> <td>12</td> <td>20cm</td> <td>038"</td> </tr> <tr> <td>14</td> <td>20cm</td> <td>038"</td> </tr> <tr> <td>11</td> <td>20cm</td> <td>038"</td> </tr> <tr> <td>6</td> <td>20cm</td> <td>038"</td> </tr> <tr> <td>9</td> <td>20cm</td> <td>038"</td> </tr> </tbody> </table>	Fr	Length	Wire Diameter	4	20cm	025"	4	20cm	035"	5	20cm	035	6	20cm	035"	6	20cm	038"	7	20cm	038"	8	20cm	038"	9	20cm	038"	10	20cm	038"	11	20cm	038"	12	20cm	038"	14	20cm	038"	11	20cm	038"	6	20cm	038"	9	20cm	038"	K15/042025	10
	Fr	Length	Wire Diameter																																																
	4	20cm	025"																																																
	4	20cm	035"																																																
	5	20cm	035																																																
	6	20cm	035"																																																
	6	20cm	038"																																																
	7	20cm	038"																																																
	8	20cm	038"																																																
	9	20cm	038"																																																
	10	20cm	038"																																																
	11	20cm	038"																																																
	12	20cm	038"																																																
	14	20cm	038"																																																
	11	20cm	038"																																																
	6	20cm	038"																																																
9	20cm	038"																																																	
		K15/042035	10																																																
		K15/052035	10																																																
		K15/062035	10																																																
		K15/062038	10																																																
		K15/072038	10																																																
		K15/082038	10																																																
		K15/092038	10																																																
		K15/102038	10																																																
		K15/112038	10																																																
		K15/122038	10																																																
		K15/142038	10																																																
		K15/17/H	10																																																
		K15/2/H	10																																																
		K15/5/H	10																																																